



APRV: A Ray of Hope Despite all of the Darkness

Authors

Dr Anant Vikram Pachisia, Dr Sandip Bhattacharya

A 60 year old male, known case of morbid obesity with obstructive sleep apnoea presented with complains of increasing breathlessness for 2 days. He was afebrile, drowsy and dyspnoeic, BMI was 44 Kg/m² (Obesity class III) with heart rate 116/min, blood pressure 170/90 mmhg, SpO₂ 64% on room air and Respiratory rate 26-30/min. On auscultation heart sounds were normal and bilateral air entry was decreased in lungs with presence of wheeze. Arterial blood gas was suggestive of type II respiratory failure with severe respiratory acidosis. Initially BIPAP was applied. However, the patient subsequently developed bradycardia. He was intubated and revived by Cardiopulmonary resuscitation as per ACLS protocol. Post resuscitation he was haemodynamically stable with GCS E₂V_TM₃.

Diagnosis of severe Acute respiratory distress syndrome (ARDS) was made based on PaO₂/FiO₂ of 78.75 as per Berlin definition.^[1] Chest x-ray showed bilateral opacities in lung parenchyma. 2D-Echo was normal. He was initially ventilated with Volume control (VC) mode with 80% FiO₂, 500ml TV, 24/min RR and PEEP of 14 cm of H₂O. He was sedated and paralysed with infusions of Inj. Fentanyl and Inj. Atracurium. However hypoxia did not improve as PaO₂ was 58 mmhg but PCO₂ decreased to 48mmhg. Peak airway pressures were 40-42 cm of H₂O.

Recruitment manoeuvre was tried with PEEP of 30 cm of H₂O, however it was not fruitful as PaO₂ was same and SpO₂ was 86% with 80% FiO₂. The mode of ventilator was changed to Pressure control (PC) with 75% FiO₂, 22/12 cm of H₂O PS/PEEP, 20/min RR. With PC mode, Peak airway pressures were similar and also PaO₂ was 55mmhg.

We considered placing the patient in prone position, as suggested by PROSEVA trial.^[2] However as the patient was morbidly obese, the decision was deferred as it would have significantly affected the nursing care. Finally the patient was ventilated by APRV mode with initial settings: 70% FiO₂, 30 cm of H₂O P high, 6 cm of H₂O P low, 5 sec T high, 0.5 sec T low. PaO₂ increased to 65 mmHg within 2 hours and SpO₂ increased to 90%. He was tracheotomised in i/v/o low GCS. By day 10 in ICU, with 55% FiO₂, PaO₂ and SpO₂ increased to 70 mmhg and 94% respectively. P high was gradually decreased to 22 cm of H₂O but P low was same. There was no need for paralysis with APRV mode and requirement for sedation was occasional. Peak pressures were 32-35 cm of H₂O.

Patient was weaned from APRV to PSV mode with PS/PEEP of 20/8 cm of H₂O and 50% FiO₂. Finally shifted to portable BIPAP machine with similar settings. We believe since muscle

relaxation was not required with APRV mode, we could wean the patient early from the ventilator. However, we ultimately lost the patient to Septic shock with multi-organ failure.

Our experience is supported by a randomized controlled trial, in which 30 mechanically ventilated trauma patients were randomly assigned to either APRV or pressure-limited ventilation. APRV was associated with shorter duration of mechanical ventilation, a shorter ICU length of stay, and use of less sedatives and paralytics.^[3]

Thus to conclude, morbidly obese patients with ARDS, in whom conventional modes of ventilation and other adjunctive therapies may be of questionable benefit can show good response to APRV mode.

Financial support and sponsorship: Nil.

Conflict of interest: There are no conflicts of interest.

References

1. The ARDS Definition Task Force. Acute Respiratory Distress Syndrome: The Berlin Definition. *JAMA*. 2012;307 (23): 2526-2533. doi:10.1001/jama.2012.5669
2. Guerin C, Reignier J, Richard JC. Prone positioning in severe acute respiratory distress syndrome. *N Engl J Med*. 2013; 368(23):2159–2168.
3. Putensen C, Zech S, Wrigge H, Zinserling J, Stuber F, Von Spiegel T, Mutz N. Long-term effects of spontaneous breathing during ventilatory support in patients with acute lung injury. *Am J Respir Crit Care Med*. 2001;164(1):43–49.